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I. TECHNICAL BASIS

This enzymatic immunoassay (ELISA) has been developed for the detection of specific antibodies to *Micobacterium avium ssp paratuberculosis* (MAP) in ruminant serum and bovine milk samples. Samples must be preincubated with a diluent with a *Micobacterium phlei* to eliminate crossreactions and then incubate in the wells coated with MAP antigen. After a first washing

step, the anti-ruminant IgG conjugate is added and after a second wash, the substrate is dispensed. The colour development indicates the presence of specific antibodies to MAP. The more concentration of antibodies it has, the more intensity of colour it develops, and therefore, the optical density (OD) will be higher.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the kit, control positive and negative serum must be tested in a systematic way.
10. **IMPORTANT!** Substrate must be handle with care, it is very sensitive to light and contamination. Remove the volume to use from the bottle and never put the remaining volume back into it.
11. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic microplate washer or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turnover of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as it is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing, shake the plate turned over on an absorbent filter paper.

V. PREPARATION OF SAMPLES:

Dilute and pre-incubate samples with M. phlei green diluent (DEMP-01). These mixes must be done in wells of empty predilution microplate or in low volume vials.

Milk:

Both whole milk and skimmed or partially skimmed milk or whey can be tested. Samples must be diluted 1/2 in diluent (60µl milk samples + 60µl diluent).

Sera:

Samples must be diluted 1/20 in diluent (6 µl serum samples + 114 µl diluent)

Diluted samples must be preincubated 30 min and then 100 µl of the mixture is transferred to the antigen plate supplied in the kit.

VI. PREPARATION OF REAGENTS

- **Washing solution:**

Dilute one part of the concentrate washing solution provided in the kit with 24 parts of distilled water or deionized water (40 ml of concentrate solution and 960 ml of water). This solution remains stable between +2°C and +8°C at least 3 months

- **Diluent of samples (M.phlei DEMP-01):**

Is ready to use. Do not dilute. Shake before to use.

- **Positive and negative controls:**

Positive and negative controls must be diluted as serum samples (6µl controls+114µl M.phlei diluent)

- **Preparation of the conjugate:**

Is ready to use. Do not dilute

VII. TEST PROCEDURE

1. All reagents must be allowed to come to room temperature before use.
2. In an empty predilution microplate dispense the following amounts of reagents:
Milk: add 60 µl of milk samples and 60 µl of M.phlei diluent.
Sera and controls: add 6 µl of controls and sera samples and 114 µl of M.phlei diluent. We recommend to run controls in duplicate wells Cover the plate and incubate **30 min ± 3 min at 22-26°C**.
3. Transfer 100 µl of mix to microplate MAP coated included in the kit. Cover and incubate the plate for **45 min ± 4min at 22-26°C**
4. Wash 3 times with washing solution following the described procedure.
5. Add 100µl of the ready to use conjugate to each well. Cover and incubate the plate for **30 min ± 3min at 22-26°C**.
6. Wash 3 times following the described procedure.
7. Add 100 µl of substrate solution, to each well. Shake softly and keep the plate for **15 min at 22-26°C** in a dark place.
8. Add 100 µl of stop solution to each well in the same order as the substrate was dispensed.
9. Read the OD of each well with a spectrophotometer at 450 nm in the following 5 minutes.

VIII. READING AND RESULT INTERPRETATION

B. Sample/positive Index Calculation (S/P %):

If you run samples in duplicate, OD values will be calculated for each sample as the arithmetic mean of both values.

To obtain S/P % of the sample:

$$S/P\% = \frac{(\text{OD Sample} - \text{OD Neg C})}{(\text{OD Pos C} - \text{OD Neg C})} \times 100$$

C. Validation of the test:

The test could be considered valid if:

1. OD_{450nm} Positive control ≥ 0.75
2. OD_{450nm} Negative control < 0.2

D. Interpretation of the results:

Milk:

- Sample with S/P % ≥ 30 will be considered as positive.
- Sample with S/P % ≤ 20 will be considered as negative.
- Sample with S/P % between both values will be considered doubtful.

Serum:

- Sample with S/P % ≥ 55 will be considered as positive.
- Sample with S/P % ≤ 45 will be considered as negative.
- Sample with S/P % between both values will be considered doubtful.

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

Eurofins-INGENASA
Avda. de La Institución Libre de Enseñanza,
39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in.....by:

IT-73840
IT-73780

ISO 14001

9191.INGE

ISO 9001:2015

9175.ING2