# INGEZIM Brucella Bovina 2.0

Prod Ref: 12.BB2.K1

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## I. TECHNICAL BASIS

Brucella species have two main antigenic components. One of those is a protein-polysaccharide complex LPS-S with endotoxic activity, responsible of the specificity of agglutination and is the major antigenic component. The second one is a polysaccharide (poli-B) with no endotoxic activity and with no relevance role in the agglutination test.

The INGEZIM BRUCELA kit is based on the indirect immunoenzymatic assay described below:

The antigen is fixed in a solid support (polystyrene plate). The antigen used is a purified extract of the

LPS of Brucella. The serum sample is added to each well. After incubation and washing we add a labelled monoclonal antibody (peroxidase conjugated) specific to bovine immunoglobulins. If the serum sample contains antibodies against the Brucella, the conjugate will bind to the plate, whereas if it does not contain specific antibodies no binding of the conjugate will be appreciated. After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled conjugate by adding the substrate that in presence of the peroxidase will develop a colorimetric reaction.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
10. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
11. Substrate must be handled with care, it is very sensible to light and contamination.

## III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

## IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF REAGENTS

- **Washing solution:**

Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water (40 ml of the concentrated solution in 960 ml of distilled or deionized water). When ready this solution remains stable between +2°C and +8°C.

- **Preparation of the conjugate (to make immediately before use):**

Dilute the needed quantity of conjugate provided in the Kit 1/100 with diluent:

- ❖ The necessary and sufficient quantity of conjugate for a complete plate is 110 µl of Conjugate in 11 ml of diluent.

- ❖ The necessary and sufficient quantity of Conjugate for an eight wells strip is 10 µl of conjugate in 1 ml of diluent.

Shake very well the solution before use. Prepare only the quantity needed for each time because the remainder volume has to be rejected.

- **Controls:**

Controls are ready to use. **DO NOT DILUTE**

## VI. PREPARATION OF SAMPLES:

### Individual serum:

Dilute 1/100 in diluent supplied. We recommend making a previous dilution 1/10.

To prepare these samples, we recommend using an empty plate (not coated one). Add 90 µl of diluent and 10 µl of each sample to each well (dilution 1/10)

### Pooled serum:

Mix equivalent volumes of the sera to be assayed (8 sera). Once homogenised, prepare a 1/25 dilution in diluent. To prepare this dilution we recommend adding 240µl of diluent and 10µl of the pool of 8 wells in an uncoated well. Homogenise adequately before use.

## VII. TEST PROCEDURE

1. All reagents (except conjugate) must be allowed to reach room temperature before use. We recommend to warm the diluent previously for 20-30 min at 37°C.

2. **Serum**

Add 90 µl of diluent in each well except the ones destined to controls. Add 10 µl of the diluted **individual samples** (diluted 1/10 as it's specified in the previous instructions) to remainder wells of the plate. Lastly, add 100 µl of the controls ready to use.

In case of **pools**, add 100 µl/well of dilution 1/25 previously prepared (point VI).

**Incubate for 1 hour at room temperature (20-25°C)**

3. Wash 5 times following the described procedure.

4. Add 100 µl of conjugate (prepared following previous instructions) to each well. Incubate the plate for **30 minutes at room temperature (20-25°C)**.

5. Wash 5 times following the described procedure.

6. Add 100 µl of substrate to each well. Keep the plate in darkness for **10 min** at room temperature.

7. Add 100 µl of stop solution to each well.

8. Read the OD of each well with a spectrophotometer at 450 nm within 5 min after the addiction of stop solution.

## VIII. READING AND RESULT INTERPRETATION

---

The reading must be done with a spectrophotometer at **450 nm**.

### A. VALIDATION OF THE TEST: SERUM

---

The test is considered valid when:

- OD value of the positive control serum is  $\geq 1.0$
- OD value of the negative control serum is  $\leq 0.2$ .

### B. INTERPRETATION OF THE RESULTS:

---

Positive and negative cut offs are calculated as follows:

#### Individual sera:

*Cut Off* = Abs 450 nm Control (+)  $\times 0.4 = 40\%$  Percentage of Positivity

#### Pools of sera:

*Cut Off* = Abs 450 nm Control (+)  $\times 0.35 = 35\%$  Percentage of Positivity

Samples with OD higher than or equal to the CUT OFF value must be considered as Positive.

Samples with OD values lower than the CUT OFF values must be considered as Negative.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tif: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: [ingenasa@ingenasa.com](mailto:ingenasa@ingenasa.com)  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in \_\_\_\_\_ by:

IT-73840  
IT-73780

ISO 14001:2015  
9191.INGE

ISO 9001:2015  
9175.ING2