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I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). The technique is briefly described below:

The recombinant antigen of *Mycobacterium bovis*, is fixed on polystyrene plates. When a serum sample contains specific antibodies to *M.bovis*, they will bind to the recombinant antigen adsorbed on the plate. After a washing step to eliminate all non-fixed material, the presence of swine antibodies can be detected using a specific peroxidase conjugate. After the

addition of the substrate, a colorimetric reaction will appear which can be measured by a spectrophotometer at 450nm.

Please note that our kit uses recombinant antigen obtained by the expression of specific and highly immunogenic *M.bovis* proteins. This method warrants the total absence of infectivity in any of the components of the kit.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions for use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix reagents or use instructions from different kits.
4. Avoid any contamination of the kit reagents.
5. Do not use components after expiry dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking while handling reagents or samples.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Systematically include a positive control and a negative serum each time the assay is run.
10. The substrate must be handled with care, as it is very sensitive to light and contamination.
11. The stop solution is a strong acid and therefore must be handled with care.

III. STORAGE OF COMPONENTS

Store all plates and reagents between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING PROCEDURE

The washing steps can be done using an automatic plate washer or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the instructions below:

- Briskly turn the plate over to empty the wells.
- Repeat the process as many times as indicated in the kit instructions.
- Prior to emptying the content of the plate in the last washing step, verify that the next reagent to be added is ready to be used. Do not let the plate dry longer than strictly necessary.
- After the last washing step tap the plate upside down on absorbent filter paper to remove any remaining washing solution.
- Throw away the content of the plate by briskly turning the plate over in order to avoid the possible exchange of the contents from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake the plate gently, avoiding cross-contamination between wells.

V. PREPARATION OF SAMPLES

- **Serum:** Serum samples must be tested at 1/200 dilution in diluent DE34-01 (i.e. 2 µl of serum + 398 µl of diluent). If the number of samples is high, it is recommended that a 1/10 dilution in an empty plate (10 µl + 90 µl serum diluent) and then transferred 5 µl of the test plate on which has previously been 95µl of diluent added (1/20 dilution)
- **Blood spots on filter paper:** See annexe

VI. PREPARATION OF REAGENTS

- **Washing solution:**
Dilute one part of the concentrated washing solution provided in the kit with 24 parts of distilled or deionized water (40 ml of concentrated solution and 960ml of water). Once prepared this solution remains stable when stored at +4°C.
- **Control serum (+) and (-):**
Both controls are ready to use. Do not dilute
- **Conjugate:**
This reagent is ready to use. Do not dilute.
- **Serum diluent:**
This reagent is ready to use. Do not dilute.

VII. TEST PROCEDURE

1. Prior to starting the test, bring all reagents to room temperature (22-25°C).
2. Add 100 µl of positive control to two wells of the plate, 100 µl of the negative control to another 2 wells and 100 µl of each sample, prepared according to previous instructions, in the remaining wells of the plate. Seal the plate and incubate for **60 minutes** at **room temperature** (18-25°C).
3. Wash 3 times following the procedure previously described.
4. Add 100 µl of the ready to use conjugate to each well. Seal the plate and incubate for **30 minutes** at **room temperature** (18-25° C).
5. Wash 5 times following the procedure previously described.
5. Add 100 µl of the substrate solution to each well. Keep the plate at **room temperature** for **10 minutes**. In order to speed up this process, it is advisable to use a multichannel pipette. Add 100 µl of the stop solution to each well. We recommend adding this reagent following the same order in which the substrate was added.
6. Read the OD of each well at 450 nm.

VIII. READING AND INTERPRETATION OF THE RESULTS:

Validation of the test:

The test is considered valid when the OD of the positive control minus the OD of the negative control is higher than 1.25 and the negative control OD is lower than 0.20.

Interpretation of the results:

The S/P must be calculated as follows:

$$S/P = \frac{(\text{Sample OD}) - (\text{Neg Control OD})}{(\text{Pos Control OD}) - (\text{Neg Control OD})}$$

Samples with S/P higher or equal than 0.25 should be considered **positive**

Samples with S/P lower than 0.25 should be considered **negative**

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Use a punch to obtain a disc about 6 mm in diameter

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
Place the disc in an empty ELISA plate or in an eppendorf tube and add 250 µl of the diluent provided in the kit

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Incubate 1 hour at 37°C or overnight (16-18 h) at 4°C without shaking

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Use the extraction to the dilution 1/2, adding first 50 µl of diluent and secondly 50 µl of the extraction on the coated plate provided in the kit.

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by:

IT-73840
IT-73780

ISO 14001:2015
9191.INGE

ISO 9001:2015
9175.ING2