INGEZIM BVD DAS

Prod Ref: 12.BVD.K2

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I. TECHNICAL BASIS

This kit is based on a double antibody sandwich enzymatic immunoassay (DAS or capture Elisa). We make a brief description of the technique below:

The plate is coated with a monoclonal antibody to the non-structural protein p80/p125 of the BVD virus. In this way when a sample is added into the wells, if it contains antigen, the antibodies will catch it. After incubation and washing, a

biotin-conjugated monoclonal antibody to p80/p125 is added.

If there is antigen captured into the well, the biotin conjugate, will bind to it. After incubation and washing you add a second conjugate (Streptavidin peroxidase conjugated) and the peroxidase action will be detected after adding the adequate substrate.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use an new tip for each serum sample.
9. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
10. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
11. Substrate must be handle with care, it is very sensible to light and contamination.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. LEUKOCYTE PREPARATION TECHNIQUE FROM WHOLE BLOOD:

a) *Materials:*

- ✓ Heparin or EDTA blood.
- ✓ 5 ml tube
- ✓ 5 ml pipette.
- ✓ Vortex
- ✓ Centrifuge
- ✓ Lysis buffer (store between +2°C and +8°C)
 - NH₄Cl : 16.6 g
 - NaH CO₃ : 2 g
 - Na₂ EDTA: 0.185 g
 - Distilled H₂O untill 1 L
 - pH 7,2-7,4
- ✓ Bicarbonate buffer 25 mM
 - NaHCO₃ : 0,21 g
 - Distilled H₂O untill 100 ml.

b) *Method*

- ✓ Centrifuge 2 ml of blood at 1000x g for 15 minutes
- ✓ Remove the supernatant (plasma) and store it frozen if needed
- ✓ Add 2 ml of lysis buffer to the cell pellet. Mix thoroughly by vortexing and leave at room temperature for 15 minutes.
- ✓ Centrifuge at 1000x g for 15 minutes.
- ✓ Discard the supernatant and add 400 µl of bicarbonate buffer to the pellet. Vortex thoroughly.
- ✓ Use 100 µl/well for testing as described in test protocol.

VI. PREPARATION OF SAMPLES:

• *Plasma, serum, whole blood samples:*

Whole blood (with heparin or EDTA), plasma or serum can be used (in this order of preference). **Before use, we recommend, although it is not mandatory to freeze and thaw the samples once to avoid false positives results.**

Samples are used without dilution (100 µl/well). Samples doubtful or from animals younger than 6 months should be tested using a leukocyte preparation (see section V). Leukocytes can be tested without freezing and thawing process.

The highest sensitivity is reached by testing leukocytes samples.

• *Ear tissue samples:*

1. Transfer the ear sample to a plate of deep wells.
2. Add 200µl of Extraction Buffer to the wells. Pressing lightly the sample with the pipette tip.
3. Incubate over night (12-18 hours) in the range of 4°C and 22°C.

• *Cell culture:*

Cell cultura samples will be assayed without dilution. In case of titration, dilute in sample diluent.

VII. PREPARATION OF REAGENTS

• *Washing solution:*

Dilute one part of the concentrate washing solution provided in the kit with 24 parts of distilled or deionized water (40 ml concentrated solution with 960 ml of distilled water). When ready this solution remains stable at +4°C

• *Controls:*

Controls are at working dilution, ready to use.

• *Diluent:*

1 PLATE FORMAT

The diluent is at working dilution, ready to use.

5 PLATES FORMAT

5x concentrated. Before using it, a 1/5 dilution must be done with distilled or desionized water (1 volume of concentrated supplied diluent with 4 volumes of distilled water)

• *Extraction buffer:*

1 PLATE FORMAT

The extraction buffer is at working dilution, ready to use.

5 PLATES FORMAT

5x concentrated. Before using it, a 1/5 dilution must be done with distilled or desionized water (1 volume of concentrated supplied diluent with 4 volumes of distilled water)

• *Preparation of the conjugates: to make immediately before use.*

Dilute the needed quantity of conjugates provided in the kit 1/100 into diluent (e.g. 10 µl of conjugate into 1 ml of diluent is enough quantity for one strip of 8 wells). For a plate add 110 µl of concentrated to 11 ml of diluent. Shake very well the solution before the use.

Prepare only the quantity needed for each time because the remainder volume has to be rejected

VIII. TEST PROCEDURE I (serum, plasma and whole blood)

- All reagents (except conjugates) must be allowed to come to room temperature before use.
- Add 100 µl/well of positive control, negative control and samples to test. We recommend the use of duplicate wells for both, samples and controls. Seal the plate and incubate **for 18 h at room temperature (18-25°C)**.
- Wash six times as specified in point IV.
Verify the cleanliness of the wells (no residual traces of blood) and if necessary, proceed with additional washes as advised in point V
- Add 100 µl/well of biotin conjugate or conjugate 1 (prepared following previous instructions) to each well. Seal the plate and **incubate for 1 hour at 37°C**
- Wash 6 times following the described procedure.
- Add 100 µl/well of peroxidase conjugate or conjugate 2 (prepared following previous instructions) to each well. Seal the plate and **incubate for 30 minutes at room temperature (18-25°C)**.
- Wash 6 times following the described procedure.
- Add 100 µl/well of TMB substrate. **Keep the plate for 5 minutes at room temperature.**
- Add 100 µl/well of stop solution.
- Read with spectrophotometer at 450nm within 5 min after the addition of stop solution.

IX. READING AND RESULT INTERPRETATION. (test procedure I)

A. Validation of the test

The test can be considered valid when:

- ⇒ The ratio, OD of the Positive control / OD of Negative control is higher than 10

B. Results interpretation:

When you are running duplicates, OD of the sample will be calculated as the arithmetic mean of OD values in both wells.

CUT OFF = OD of Positive Control x 0.1

- Samples with an OD higher than cut off value plus 15% must be considered as **POSITIVE** to BVDV antigen
- Samples with OD lower than cut off value minus 15% must be considered as **NEGATIVE** to BVDV antigen
- Samples with OD values between both values are considered as **DOUBTFUL**. We recommend re-testing these animals one more time or applying a different technique to check this sample.

C. Example:

OD Control (+) = 2.05

OD Control (-) = 0.07

Cut off = $2.05 \times 0.1 = 0.205$

- Samples Positives (OD higher than)
 $0.205 + (0.15 \times 0.205) = 0.205 + 0.030 = 0.235$
- Samples Negative (OD lower than)
 $0.205 - (0.15 \times 0.205) = 0.205 - 0.030 = 0.174$
- Samples doubtful
OD between 0.174 and 0.235

Any positive sample should be considered as sample from a persistently infected animal. However, you need two positive results of paired samples of the animal (interval of 3-4 weeks between them) in order to confirm the diagnosis. In addition, if the animal is antibody-negative to BVDV you must definitively consider it as PI.

X. TEST PROCEDURE II (ear tissue and cell culture)

1. All reagents (except conjugates) must be allowed to come to room temperature before use.
2. Add 100 µl/well of positive control, negative control and samples to test. We recommend the use of duplicate wells for both, samples and controls. Seal the plate. We recommended to incubate **for 2,5 h at 37°C**. Nevertheless, there is the possibility of incubate over night between +2°C and +8°C).
3. Wash six times as specified in point IV.
4. **Verify the cleanliness of the wells (no residual traces of blood) and if necessary, proceed with additional washes as advised in point V**
5. Add 100 µl/well of biotin conjugate or conjugate 1 (prepared following previous instructions) to each well. Seal the plate and **incubate for 1 hour at 37°C**.
6. Wash 6 times following the described procedure.
7. Add 100 ul/well of peroxidase conjugate or conjugate 2 (prepared following previous instructions) to each well. Seal the plate and **incubate for 30 minutes at room temperature (18-25°C)**.
8. Wash 6 times following the described procedure.
9. Add 100 µl/well of TMB substrate. **Keep the plate for 5 minutes at room temperature**.
10. Add 100 µl/well of stop solution.
11. Read with spectrophotometer at 450nm within 5 min after the addition of stop solution.

XI. READING AND RESULT INTERPRETATION. (test procedure II)

A. Validation of the test

The test can be considered valid when:

- ⇒ The ratio, OD of the Positive control / OD of Negative control is higher than 10

B. Results interpretation:

When you are running duplicates, OD of the sample will be calculated as the arithmetic mean of OD values in both wells.

CUT OFF = 10% of OD of Positive Control

- Samples with an OD higher than or equal to cut off value must be considered **POSITIVE** to BVDV antigen
- Samples with OD lower than cut off value must be considered as **NEGATIVE** to BVDV antigen

C. Example:

OD Control (+) = 1.65

OD Control (-) = 0.07

Cut off = $1,65 \times 0.1 = 0.165$

- Samples Positives: OD higher than or equal to 0.165
- Samples Negative: OD lower than 0.165

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by: