INGEZIM BLV COMPAC 2.0

Prod Ref: 12.BLV.K3

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I. TECHNICAL BASIS

This Kit is based on a blocking enzymatic immunoassay (Blocking Elisa) based on the use of two monoclonal antibodies against the viral gp51. This test allows the detection of antibodies against the BLV in bovine serum or milks (individual or pooled). We make a brief description of the technique below:

The solid phase are plates coated with BLV gp51, bound to the plate through two specific monoclonal antibodies against the gp51 virus protein. After adding the sample to the well, if it contains specific antibodies against the virus, they will bind to the antigen absorbed on plate while if the sample does not contain specific antibodies they will not bind to the antigen.

If we add a specific monoclonal antibody against the viral antigen coated to the plate (conjugated with peroxidase), it will compete with the antibodies of the serum. If the serum samples contains specific antibodies, they will not permit the binding of the labelled Mab to the antigen whereas if it does not contain specific antibodies the Mab will bind to the antigen on the plate. After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled Mab by adding the substrate (TMB) that in presence of the peroxidase will develop a colorimetric reaction.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
10. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
11. Substrate must be handled with care, it is very sensible to light and contamination.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well..
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF REAGENTS

- **Washing solution:**

Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water(40 ml of concentrate solution and 960 ml of water). Once this solution ready, it remains stable between +2°C and +8°C for 3 months maximum as long as adequate conditions will be maintained.

- **Positive and Negative controls:**

Controls should be handled like samples, adding 50 µl/well.

- **Preparation of the conjugate (to make immediately before use):**

Dilute the needed quantity of conjugate

provided in the Kit 1/100 with diluent::

⇒ The necessary and sufficient quantity of conjugate for a complete plate is 110 µl of Conjugate in 11 ml of diluent.

⇒ The necessary and sufficient quantity of Conjugate for an eight wells strip is 10 µl of conjugate in 1 ml of diluent.

Shake very well the solution before use.

Prepare only the quantity needed for each time because the remainder volume has to be rejected.

VI. PREPARATION OF SAMPLES:

- **Use of serum:**

Samples can be tested individually or pooled (up to 10 serum samples).

- **Use of milk:**

Samples can be tested only individually. It is recommended to skim the milk samples before testing.

VII. TEST PROCEDURE

1. All reagents (except conjugate) must be allowed to come to room temperature before use.
2. Add 50 µl of diluent DE04-01 to each well of the plate to be used. Add 50 µl of Positive control, Negative control and samples to the wells of the plate. In this way you are assaying the samples at a 1/2 dilution. Shake the plate carefully for the better homogenization of the dilution. We recommend running control sera in duplicate using 2 wells for each one. In order to avoid border effect detected when some kind of incubators are used, it is recommendable not use the external files to add the controls). **Incubate the plate 1-hour at 37°C.**
3. Wash 3 times following the described procedure.
4. Add 100 µl of conjugate (prepared following previous instructions) to each well. **Incubate the plate for 30 minutes at +37°C.**
5. Wash 4 times following the described procedure.
6. Add 100 µl of substrate solution, to each well. Keep the plate for **10 min at room temperature**, in a dark place.
7. Add 100 µl of stop solution to each well following the same order in which the substrate was added.
8. Read the OD of each well with a spectrophotometer at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm**.

* If the OD value of the Negative Control Serum is more than 5 times the OD value of the Positive control Serum:

A. Validation of the Test:

The test could be considered valid:

$$\frac{\text{OD (Negative Control Serum)}}{\text{OD (Positive Control Serum)}} > 5$$

- * The OD of the Negative Control Serum must be > than 1.

B. Interpretation of the Results:

1. **Cut Off calculation:**

$$\sum \text{Negative Cut Off} = \text{NC} - [(\text{NC} - \text{PC}) \times 0.4]$$

$$\sum \text{Positive Cut Off} = \text{NC} - [(\text{NC} - \text{PC}) \times 0.5]$$

Where: NC = Mean of OD values of the Negative control serum.
PC = Mean of OD values of Positive Control Serum.

2. **Results Interpretation:**

When you are running duplicate samples, OD values will be calculate for each sample as the arithmetic mean of both values.

⇒ All samples with OD higher than or equal to negative CUT OFF value must be considered as Negatives.

⇒ All samples with OD values lower than or equal to Positive CUT OFF values will be considered as Positives

Samples with OD values between both cut-offs must be considered as doubtful. For these samples a new assay is recommended after 3 weeks.

EXAMPLE:

Means of OD

* Positive Control (PC)	0.154
* Negative Control (NC)	1.750
* Sample 1	0.170
* Sample 2	1.684

1. **Test Validation:**

$$\text{NC/PC} = 11.3 > 5$$

$$\text{NC} = 1.750 > 1$$

2. **Cut off calculation:**

$$\text{CUT OFF (-)} = 1.750 - [(1.750 - 0.154) \times 0.4] = 1.111$$

$$\text{CUT OFF (+)} = 1.750 - [(1.750 - 0.154) \times 0.5] = 0.952$$

3. **Results Interpretation:**

Sample 1 » Positive to BLV antibodies
Sample 2 » Negative to BLV antibodies

For the interpretation of results a Microsoft excel file is available in our web.

Developed and manufactured in Spain by:

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INMUNOLOGIA Y GENETICA APLICADA, S.A.